

PROST "save"

Medical device for the prevention of secondary effects of localized compression in the pelvic area



CHAIRS

Frequent Usage

- » Interventions to the urogenital apparatus in men and in women
- » Micro fracture of coccyx
- » Hemorrhoids
- » Prevention of BPH (Benign Prostatic Hypertrophy)
- » Postpartum
- » Sedentariness

Prevention & Convalescence



www.gardenbalance.it



what is?

A special chair...or rather a medical device that can be used for the prevention of secondary effects of localized compression in the pelvic area. Unique in its kind, **Prost "Save"** is designed to provide relief and well-being in cases of post-surgical convalescence, chronic diseases, trauma, injury, inflammation, sedentary, etc., affecting the pelvic area. **Prost "Save"** more efficiently replaces the numerous removable pillows whose positioning is less comfortable and not always adequate.



Also particularly suitable in the presence of catheters and drains.

Do you want to change the color of the upholstery?

Frame or click the QR codes and find out how to customize your device by browsing our color cards*.


*If no preference is expressed, the device will be supplied with the available color.






How It Acts

Direct contact of the organs of the pelvic area with a flat, albeit soft surface (common chair - image 1), generates static compressions that, in the long term, can be the cause or with cause of many severe pathologies. As a mere example, in men, crushing of the prostate (an outcome of sedentariness), can cause sometimes degenerative inflammation.



The distribution of body mass that does not affect the pelvic area (image 2) is indicated for the prevention of any possible pathology affecting the urogenital system and rectum.





When to use it...



... POST SURGERY

is useful for avoiding compression in the pelvic area after uro-genital and rectal surgeries, regardless of gender. The use of this medical device can promote, and accelerate, healing.

The special geometry of the Prost "Save" chair is compatible with the presence of catheters and drains.

... FRACTURE OF THE COCCYX

to decompress the sacral area ... to accelerate the consolidation of the bone injury and achieve an antalgic effect.



... PRESSURE SORES

for the prevention of inflammation and pressure sores in the sacral area.

... HEMORRHOIDS

for prevention and treatment of irritation and hemorrhoids. Soothing in contact dermatitis in the pelvic area.





... PRE AND POSTPARTUM



during pregnancy as weight gain compresses the pelvic floor causing annoying consequences...even after delivery in the presence of lacerations or episiotomy.

... SEDENTARINESS

one is forced to sit for long periods of time, the use of the prost "Save" chair proves functional for the prevention of the related complications (prostatitis, bladder inflammation, skin lesions, etc.).





Benefits

- because the cushion is integral to the chair, accommodation is extremely practical and safe;
- padding with multiple overlapping layers of differentiated density polyurethane foam and memory foam on the surface;
- upholstery available in Similpelle (standard equipment) and Silica® (optional component - code AC87) - upholstery color of your choice;
- surfaces that can be sanitized with commonly used detergents;
- medical device designed to withstand heavy use;
- usable by anyone, anywhere;
- low weight and large capacity;
- small footprint;
- chrome-plated steel frame;
- non-slip feet;
- ergonomic finishes.





Data sheet

CND and EMDN	V080299
UDI-DI di Base	805771740DITRGR
No. progressive/R	2564199/R
ID Product	GB0428.SP
Intended Use	Devices intended for patient accommodation during hospitalization and medical examinations and procedures in general
Max load	130kg

Dimensions

Backrest Height	85 cm
Seat Height	55 cm
Width	52 cm
Depth	55 cm

For different sizes and capacities, request a quote.



PROST "save"

Medical device for the prevention of secondary effects of localized compression in the pelvic area

www.gardhenbalance.it



Gardhen Balance S.r.l

Via G. Luraghi c/o Consorzio Il Sole Lotto S
80038 Pomigliano d'Arco (NA) - Italy
Tel. +39.081.8692160 pbx
Fax +39.081.8692460
Email: info@gardhenbalance.it

Images reproduced in this sheet are for reference only and dimensional values are nominal. It is the manufacturer's right to make changes to the designs without prior notice. Refer to the optional components catalog and/or commercial proposals for identification of standard equipment and selection of any additional components. In case of discrepant data, the contents of the commercial proposals shall prevail. For regulatory references, please refer to the user manual and declaration of conformity.
Dep. EN_382 rev. 15/04/2024